


Guiding Catheter

Used to place biliary stent endoscopically.

Technical Information

- Size - 5 fr. - 6 fr. - 7 fr.
- Working Length - 45 cm, 220 cm & 320 cm
- Availability - Sterile - Single use only.
- Individually packed.



Ordering Information-Mention product code x quantity required

Reference	Product Name	Size	Working Length	Compatible Pusher
SGC-645	GUIDING CATHETER	6FR	45 CM	-
SGC-5320	GUIDING CATHETER	5FR	320 CM	8.5 FR
SGC-5220	GUIDING CATHETER	5FR	220 CM	8.5 FR
SGC-6320	GUIDING CATHETER	6FR	320 CM	10 FR
SGC-6220	GUIDING CATHETER	6FR	220 CM	10 FR
SGC-7320	GUIDING CATHETER	7FR	320 CM	11.5 FR
SGC-7220	GUIDING CATHETER	7FR	220 CM	11.5 FR